



# ACETIAM PRI

DICOM Conformance Statement

# ACETIAM PRI

## DICOM Conformance Statement

Product Version: [5.70](#)

Documentation Update: [v18 - December 2019](#)

## Table of Contents

|   |    |
|---|----|
| ■ 1. Introduction .....   | 1  |
| 1.1 Scope and field of application .....  | 1  |
| 1.2 Acronyms and Abbreviations .....  | 1  |
| 1.3 References.....   | 2  |
| 1.4 Intended audience.....  | 2  |
| 1.5 Warning to the Reader.....  | 2  |
| ■ 2. Document history .....   | 3  |
| ■ 3. Implementation Model .....   | 3  |
| 3.1 Application Data Flow Diagram .....   | 3  |
| ■ Figure 3.1: ACETIAM PRI Implementation Model.....                                   | 4  |
| 3.2 Functional Definitions of AE's .....  | 4  |
| 3.2.1 Verification Service as SCP .....   | 4  |
| 3.2.2 Print Services as SCP .....   | 4  |
| 3.2.3 Storage Services as SCP.....  | 4  |
| 3.2.4 Verification Service as SCU.....  | 5  |
| 3.2.5 Storage Service as SCU.....   | 5  |
| 3.2.6 Modality Worklist as SCU.....   | 5  |
| 3.2.7 Query/Retrieve as SCU .....   | 5  |
| 3.3 Sequencing of Real-World Activities .....   | 5  |
| ■ 4. Application Entity Specifications.....   | 6  |
| ■ Table 4a: Supported SOP Classes for ACETIAM PRI – VERIFICATION as SCU and SCP ..... | 6  |
| ■ Table 4b: Supported SOP Classes for ACETIAM PRI – PRINT as SCP .....                | 6  |
| ■ Table 4c: Supported SOP Classes for ACETIAM PRI – STORAGE as SCP.....               | 6  |
| ■ Table 4d: Supported SOP Classes for ACETIAM PRI – STORAGE as SCU .....              | 7  |
| ■ Table 4e: Supported SOP Classes for ACETIAM PRI – MODALITY WORKLIST as SCU .....    | 7  |
| ■ Table 4f: Supported SOP Classes for ACETIAM PRI – QUERY/RETRIEVE as SCU.....        | 7  |
| 4.1 Association Establishment Policies .....  | 7  |
| 4.1.1 General.....  | 7  |
| 4.1.2 Number of Associations.....   | 8  |
| 4.1.3 Asynchronous Nature.....  | 8  |
| 4.1.4 Implementation Identifying Information.....                                     | 9  |
| ■ Table 4.1.4 : Application Identification Information for ACETIAM PRI .....          | 9  |
| 4.2 Association Initiation by Real-World Activity .....                               | 9  |
| 4.3 Association Acceptance Policy .....   | 9  |
| 4.3.1 Real World Activity .....   | 9  |
| ■ Table 4.3.1.2: Presentation Contexts.....   | 10 |
| ■ 5. Communication Profiles.....  | 19 |
| 5.1 Supported Communications Stacks .....   | 19 |
| 5.2 TCP/IP Stack.....   | 19 |
| 5.3 Physical Media Support .....  | 19 |

- 6. Extensions/Specialization/Privatization ..... 19
- 7. Configuration ..... 19
  - 7.1 AE Title/Presentation Address Mapping ..... 19
  - 7.2 Configurable Parameters ..... 19
    - 7.2.1 Standard Configuration ..... 19
    - 7.2.2 Optional Configuration ..... 20
- 8. Support of Extended Character Sets ..... 20

# 1. Introduction

## 1.1 Scope and field of application

This document describes ACETIAM PRI conformance to the DICOM 3.0 standard.

This conformance statement describes the conformance specifics of ACETIAM PRI software configurations supporting various Windows compatible devices.

It contains a short description of application involved and provides technical information about data exchange capabilities of the equipment. The main elements describing these capabilities are the supported DICOM Service Object Pair (SOP) Classes, Roles, Information Object Definitions (IOD) and Transfer Syntaxes.

It should be read in conjunction with the DICOM standard and its addenda.

This statement is conformant with the recommended format as described in PS 3.2 of the DICOM standard.

ACETIAM PRI acts as a SCP for the following Meta SOP Classes:

- Basic Grayscale Print Management
- Basic Color Print Management

ACETIAM PRI also acts as a SCP for the following SOP Classes:

- Verification SOP Class
- Presentation LUT SOP Class

ACETIAM PRI may also be configured to act as a Storage SCP to either perform automatic Store to Print operations, or to fill Print Job information with Patient demographics information found in the incoming stored objects.

## 1.2 Acronyms and Abbreviations

The following acronyms and abbreviations are used in this document

- ACR            American college of Radiology
- ANSI           American National Standards Institute
- DICOM        Digital Imaging and Communication in Medicine
- DIMSE        DICOM Message Service Element
- DIMSE-C     DICOM Message Service Element-Composite
- DIMSE-N     DICOM Message Service Element-Normalized
- NEMA         National Electrical Manufacturers Association
- PDU           Protocol Data Unit
- SCP            Service Class Provider
- SCU            Service Class User
- SOP            Service Object Pair
- TCP/IP        Transmission Control Protocol/Internet Protocol
- UID            Unique Identifier

### **1.3 References**

[DICOM]

The Digital Imaging and Communications in Medicine (DICOM) standard:

NEMA PS 3.1 - 3.20 (2011), Digital Imaging and Communications in Medicine (DICOM) Set

National Electrical Manufacturers Association (NEMA) - Publication Sales - 1300 N. 17th Street, Suite 1847 - Rosslyn, Va. 22209, United States of America.

### **1.4 Intended audience**

This Conformance Statement is intended for:

- Potential users
- System integrators of medical equipment
- Software designers implementing DICOM interfaces

### **1.5 Warning to the Reader**

It is assumed that the reader is familiar with the DICOM standard.

If another device matches this Conformance Statement based on the comparison with its own Conformance Statement, there is a chance, but no guarantee that they interoperate. DICOM only deals with communication; it is not a standard which specifies what is needed for certain applications to run on a device.

## 2. Document history

| Version | Changes     | Authors                       | Date        |
|---------|-------------|-------------------------------|-------------|
| 1       | Creation    | Fabien BOISSE                 | 16 Jul 2001 |
| 2       | Update      | Gilles MEVEL                  | 12 Feb 2002 |
| 3       | Update      | Gilles MEVEL / Jérôme GUIGNOT | 31 Mar 2003 |
| 4       | Update      | Gilles MEVEL / Jérôme GUIGNOT | 01 Jun 2007 |
| 5       | Update      | Gilles MEVEL                  | 10 Mar 2008 |
| 6       | Update      | Gilles MEVEL                  | 20 Jun 2008 |
| 7       | Update 4.22 | Gilles MEVEL                  | 11 Dec 2008 |
| 8       | Update 4.40 | Gilles MEVEL                  | 28 Sep 2009 |
| 9       | Update 4.60 | Gilles MEVEL                  | 19 Nov 2010 |
| 10      | Update 4.62 | Gilles MEVEL                  | 13 Nov 2010 |
| 11      | Update 4.64 | Gilles MEVEL                  | 20 Sep 2011 |
| 12      | Update 4.66 | Gilles MEVEL                  | 04 Mar 2013 |
| 13      | Update 5.00 | Gilles MEVEL                  | 21 Jan 2014 |
| 14      | Update 5.20 | Gilles MEVEL                  | 13 Mar 2015 |
| 15      | Update 5.30 | Gilles MEVEL                  | 11 Jan 2017 |
| 16      | Update 5.40 | Gilles MEVEL                  | 29 Jun 2017 |
| 17      | Update 5.60 | Gilles MEVEL                  | 08 Nov 2018 |
| 18      | Update 5.70 | Gilles MEVEL                  | 19 Dec 2019 |


## 3. Implementation Model

### 3.1 Application Data Flow Diagram

ACETIAM PRI is a Windows application that allows transforming your PC into a DICOM 3.0 Print Server.

ACETIAM PRI application will appear in the system tray of Windows task bar and run continuously until stopped.

**Figure 3.1: ACETIAM PRI Implementation Model**



### 3.2 Functional Definitions of AE's

ACETIAM PRI can be called with different AE Titles. This only allows ACETIAM PRI to drive several printer devices, or to have different behaviors according to the (calling AE Title, called AE Title) pair value.

#### 3.2.1 Verification Service as SCP

ACETIAM PRI waits for another application to connect at the presentation address configured for its Application Entity Title. When another application connects, ACETIAM PRI expects it to be a DICOM application. ACETIAM PRI will accept associations with Presentation Contexts for Verification Service Class.

#### 3.2.2 Print Services as SCP

Once started, ACETIAM PRI waits for association requests, and, for each accepted request, processes on the association the received print commands compatible with the SOP Classes it supports. Associations are released either on Print SCU request or when an error condition occurs that leads to an association abort.

#### 3.2.3 Storage Services as SCP

ACETIAM PRI may be configured to accept store requests, according to the calling AE Title and the called AE Title or both. When storage requests are accepted, images are received via C-STORE



commands, and stored for later automatic printing. Associations are released either on Print SCU request or when an error condition occurs that leads to an association abort.

### ***3.2.4 Verification Service as SCU***

When user performs software configuration, ACETIAM PRI Verification requests to DICOM peers.

Also, when started, ACETIAM PRI may be defined to issue DICOM Verification requests to itself, for self-testing purposes.

### ***3.2.5 Storage Service as SCU***

ACETIAM PRI may be configured to save all its production as DICOM Secondary Capture objects of as DICOM PDF files. Those DICOM objects are then sent automatically to a remote DICOM server via DICOM Store operations.

### ***3.2.6 Modality Worklist as SCU***

In some situations, when an Accession Number or a Patient ID is decoded using OCR techniques, ACETIAM PRI may be configured to perform queries to a RIS, to get study/patient related information about scheduled studies concerning the Accession Number or the Patient ID. Queries are made by sending C-FIND commands using the Modality Worklist service.

### ***3.2.7 Query/Retrieve as SCU***

In some situations, when an Accession Number or a Patient ID is decoded using OCR techniques, ACETIAM PRI may be configured to perform queries to a PACS, to get study/patient related information about performed studies concerning the Accession Number or the Patient ID. Queries are made by sending C-FIND commands using the Query and Retrieve service.

## **3.3 Sequencing of Real-World Activities**

Not Applicable.

## 4. Application Entity Specifications

ACETIAM PRI exists as a single Application Entity (AE).

In Advanced Edition, ACETIAM PRI determines its capabilities according to the calling AE Title and the called AE Title specified in the incoming association request. This section provides information on the maximum capabilities of ACETIAM PRI AE.

ACETIAM PRI provides Standard Conformance to the following DICOM SOP Classes as an SCP specified in Tables 4.1a to 4.1f:

**Table 4a: Supported SOP Classes for ACETIAM PRI – VERIFICATION as SCU and SCP**

| SOP Class Name         | SOP Class UID     |
|------------------------|-------------------|
| Verification SOP Class | 1.2.840.10008.1.1 |

**Table 4b: Supported SOP Classes for ACETIAM PRI – PRINT as SCP**

| SOP Class Name                                  | SOP Class UID           |
|---|-------------------------|
| Basic Grayscale Print Management Meta SOP Class | 1.2.840.10008.5.1.1.9   |
| Basic Color Print Management Meta SOP Class     | 1.2.840.10008.5.1.1.18  |
| Basic Film Session                              | 1.2.840.10008.5.1.1.1   |
| Basic Film Box                                  | 1.2.840.10008.5.1.1.2   |
| Basic Grayscale Image Box                       | 1.2.840.10008.5.1.1.4   |
| Basic Color Image Box                           | 1.2.840.10008.5.1.1.4.1 |
| Printer   | 1.2.840.10008.5.1.1.16  |
| Presentation LUT SOP Class                      | 1.2.840.10008.5.1.1.23  |

**Table 4c: Supported SOP Classes for ACETIAM PRI – STORAGE as SCP**

| SOP Class Name  | SOP Class UID                 |
|---|-------------------------------|
| Computed Radiography Image Storage                      | 1.2.840.10008.5.1.4.1.1.1     |
| Digital XRay Image Storage For Presentation             | 1.2.840.10008.5.1.4.1.1.1.1   |
| Digital XRay Image Storage For Processing               | 1.2.840.10008.5.1.4.1.1.1.1.1 |
| Digital Mammography XRay Image Storage For Presentation | 1.2.840.10008.5.1.4.1.1.1.2   |
| Digital Mammography XRay Image Storage For Processing   | 1.2.840.10008.5.1.4.1.1.1.2.1 |
| Digital Intra Oral XRay Image Storage For Presentation  | 1.2.840.10008.5.1.4.1.1.1.3   |
| Digital Intra Oral XRay Image Storage For Processing    | 1.2.840.10008.5.1.4.1.1.1.3.1 |
| CT Image Storage  | 1.2.840.10008.5.1.4.1.1.2     |
| <i>US Multiframe Image Storage (RET)</i>                | 1.2.840.10008.5.1.4.1.1.3     |
| US Multiframe Image Storage                             | 1.2.840.10008.5.1.4.1.1.3.1   |
| MR Image Storage  | 1.2.840.10008.5.1.4.1.1.4     |
| <i>US Image Storage (RET)</i>                           | 1.2.840.10008.5.1.4.1.1.6     |
| US Image Storage  | 1.2.840.10008.5.1.4.1.1.6.1   |
| Secondary Capture Image Storage                         | 1.2.840.10008.5.1.4.1.1.7     |
| Multiframe Secondary Capture Single Bit Image Storage   | 1.2.840.10008.5.1.4.1.1.7.1   |
| Multiframe Secondary Capture Byte Image Storage         | 1.2.840.10008.5.1.4.1.1.7.2   |

| SOP Class Name  | SOP Class UID                    |
|---|----------------------------------|
| Multiframe Secondary Capture Word Image Storage       | 1.2.840.10008.5.1.4.1.1.7.3      |
| Multiframe Secondary Capture True Color Image Storage | 1.2.840.10008.5.1.4.1.1.7.4      |
| XRay Angiographic Image Storage                       | 1.2.840.10008.5.1.4.1.1.12.1     |
| XRay Fluoroscopy Image Storage                        | 1.2.840.10008.5.1.4.1.1.12.2     |
| VL Endoscopic Image Storage                           | 1.2.840.10008.5.1.4.1.1.77.1.1   |
| VL Microscopic Image Storage                          | 1.2.840.10008.5.1.4.1.1.77.1.2   |
| VL Slide Coordinates Microscopic Image Storage        | 1.2.840.10008.5.1.4.1.1.77.1.3   |
| VL Photographic Image Storage                         | 1.2.840.10008.5.1.4.1.1.77.1.4   |
| Ophthalmic Photography 8Bit Image Storage             | 1.2.840.10008.5.1.4.1.1.77.1.5.1 |
| Ophthalmic Photography 16Bit Image Storage            | 1.2.840.10008.5.1.4.1.1.77.1.5.2 |
| Basic Text SR Storage                                 | 1.2.840.10008.5.1.4.1.1.88.11    |
| Enhanced Text SR Storage                              | 1.2.840.10008.5.1.4.1.1.88.22    |
| Comprehensive Text SR Storage                         | 1.2.840.10008.5.1.4.1.1.88.33    |
| XRay Radiation Dose SR Storage                        | 1.2.840.10008.5.1.4.1.1.88.67    |
| <i>Fuji CR Image Storage (private)</i>                | 1.2.392.200036.9125.1.1.2        |

**Table 4d: Supported SOP Classes for ACETIAM PRI – STORAGE as SCU**

| SOP Class Name                  | SOP Class UID                 |
|---------------------------------|-------------------------------|
| Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7     |
| Encapsulated PDF Storage        | 1.2.840.10008.5.1.4.1.1.104.1 |

**Table 4e: Supported SOP Classes for ACETIAM PRI – MODALITY WORKLIST as SCU**

| SOP Class Name                             | SOP Class UID          |
|--|------------------------|
| Modality Worklist Information Model - FIND | 1.2.840.10008.5.1.4.31 |

**Table 4f: Supported SOP Classes for ACETIAM PRI – QUERY/RETRIEVE as SCU**

| SOP Class Name                                     | SOP Class UID               |
|--|-----------------------------|
| Study Root Query/Retrieve Information Model – FIND | 1.2.840.10008.5.1.4.1.2.2.1 |

## 4.1 Association Establishment Policies

### 4.1.1 General

Minimum PDU size accepted: 4096 bytes

Maximum PDU size accepted: 131072 bytes

### ***4.1.2 Number of Associations***

The number of simultaneous associations supported may be defined in ACETIAM PRI configuration panel. Default is 25

### ***4.1.3 Asynchronous Nature***

ACETIAM PRI does not support asynchronous operations.

#### **4.1.4 Implementation Identifying Information**

ACETIAM PRI will respond with the following implementation identifying parameters:

**Table 4.1.4 : Application Identification Information for ACETIAM PRI**

| <b>Name</b>                 | <b>SOP Class UID</b>   |
|-----------------------------|------------------------|
| Implementation Class UID    | 1.2.250.1.59.3.0.3.5.3 |
| Application Context Name    | 1.2.840.10008.3.1.1.1  |
| Implementation Version Name | ETIAM_DCMTK_353        |

#### **4.2 Association Initiation by Real-World Activity**

For each SCU operation, ACETIAM PRI will open an association with the remote equipment, send a command, wait for its responses and close the association.

#### **4.3 Association Acceptance Policy**

According to its configuration, ACETIAM PRI may accept or not associations from remote equipment. This may depend on the Calling AETitle, the called AETitle or both.

##### **4.3.1 Real World Activity**

###### **4.3.1.1 Associated Real-World Activity**

The application entity waits for incoming associations. No operator action is required to receive ACETIAM PRI jobs or verification requests.

### 4.3.1.2 Presentation Context Table

ACETIAM PRI supports the following Presentation Contexts:

**Table 4.3.1.2: Presentation Contexts**

| Presentation Context Table  |                             |  |  |      |                      |
|---|-----------------------------|--|--|------|----------------------|
| Abstract Syntax   |                             | Transfer Syntax  |  | Role | Extended Negotiation |
| Name  | UID                         | Name   | UID                                      |      |                      |
| <i>Verification SOP Class</i>                                       | 1.2.840.10008.5.1.1         | DICOM Implicit VR Little Endian<br>DICOM Explicit VR Little Endian | 1.2.840.10008.1.2<br>1.2.840.10008.1.2.1 | SCP  | None                 |
| <i>Basic Grayscale Print Management Meta SOP Class</i>              | 1.2.840.10008.5.1.1.9       | DICOM Implicit VR Little Endian<br>DICOM Explicit VR Little Endian | 1.2.840.10008.1.2<br>1.2.840.10008.1.2.1 | SCP  | None                 |
| <i>Basic Color Print Management Meta SOP Class</i>                  | 1.2.840.10008.5.1.1.18      | DICOM Implicit VR Little Endian<br>DICOM Explicit VR Little Endian | 1.2.840.10008.1.2<br>1.2.840.10008.1.2.1 | SCP  | None                 |
| <i>Presentation LUT SOP Class</i>                                   | 1.2.840.10008.5.1.1.23      | DICOM Implicit VR Little Endian<br>DICOM Explicit VR Little Endian | 1.2.840.10008.1.2<br>1.2.840.10008.1.2.1 | SCP  | None                 |
| <i>Storage SOP Class</i>  | See table 4c                | DICOM Implicit VR Little Endian<br>DICOM Explicit VR Little Endian | 1.2.840.10008.1.2<br>1.2.840.10008.1.2.1 | SCP  | None                 |
| <i>Storage SOP Class</i>  | See table 4d                | DICOM Implicit VR Little Endian<br>DICOM Explicit VR Little Endian | 1.2.840.10008.1.2<br>1.2.840.10008.1.2.1 | SCU  | None                 |
| <i>Modality Worklist Information Model - FIND SOP Class</i>         | 1.2.840.10008.5.1.4.31      | DICOM Implicit VR Little Endian<br>DICOM Explicit VR Little Endian | 1.2.840.10008.1.2<br>1.2.840.10008.1.2.1 | SCU  | None                 |
| <i>Study Root Query/Retrieve Information Model – FIND SOP Class</i> | 1.2.840.10008.5.1.4.1.2.2.1 | DICOM Implicit VR Little Endian<br>DICOM Explicit VR Little Endian | 1.2.840.10008.1.2<br>1.2.840.10008.1.2.1 | SCU  | None                 |

### 4.3.1.3 SOP Specific Conformance for Printer SOP Class

The Printer SOP Class is used to monitor the status of the printer.

The following DIMSE services are supported:

N-EVENT-REPORT

N-GET

➔ N-EVENT-REPORT is used to report the changes of the printer status in an asynchronous way.

The attributes of N-EVENT-REPORT are shown in the following table:

| Tag         | Name                | Comment |
|-------------|---------------------|---------|
| (2110,0010) | Printer Status      |         |
| (2110,0020) | Printer Status Info |         |

The printer status will be returned as a combination of the Printer Status attribute (2110,0010) and the Printer Status Info attribute (2110,0020) of the Printer SOP Class.

| Printer Status | Printer Status Info | Meaning  |
|----------------|---------------------|--|
| NORMAL         |                     | Printer OK   |
| WARNING        |                     | General printer warning, no specific information is available. |
| WARNING        | SUPPLY EMPTY        | Printer dispenser is empty.                                    |
| WARNING        | PRINTER OFFLINE     | Printer is Offline   |
| WARNING        | FILM JAM            | No more paper or film in the Printer                           |
| FAILURE        |                     | General failure of the printer.                                |

➔ N-GET retrieves an instance of the Printer SOP Class.

The attributes of N-GET are shown in the following table:

| Tag         | Name                    | Comment   |
|-------------|-------------------------|---|
| (2110,0010) | Printer Status          | Always returns NORMAL                           |
| (2110,0020) | Printer Status Info     | Always returns NORMAL                           |
| (2110,0020) | Printer Name            | Returns the value configured in file DcmPri.cfg |
| (2110,0020) | Manufacturer            | Returns the name configured in file DcmPri.cfg  |
| (2110,0020) | Manufacturer Model Name | Returns the name configured in file DcmPri.cfg  |
| (2110,0020) | Device Serial Number    | Returns the name configured in file DcmPri.cfg  |
| (2110,0020) | Software versions       | Returns the name configured in file DcmPri.cfg  |
| (2110,0020) | Date Last Calibration   | Returns the name configured in file DcmPri.cfg  |
| (2110,0020) | Time Last Calibration   | Returns the name configured in file DcmPri.cfg  |

The SCP will return one of the following status codes for N-GET :

| Code   | Status  | Meaning           |
|--------|---------|-------------------|
| 0x0000 | Success | Success           |
| 0x0105 | Failure | No such attribute |

#### 4.3.1.4 SOP Specific Conformance for Basic Film Session SOP Class

The following DIMSE services are supported:

N-CREATE

N-SET

N-DELETE

N-ACTION

➔ N-CREATE is sent by the SCU AE to create a Basic Film Session SOP instance, when an association has been established. If N-CREATE operation fails, an error message will be returned by the SCP AE. The N-CREATE causes the Basic Film Session to be created and its attributes initialized.

The Basic Film Session SOP instances shall be created before the Film Box SOP Instances are created.

ACETIAM PRI provides the following support for the attributes contained in the N-CREATE DIMSE Service of the Basic Film Session SOP Class:

| Tag         | Name               | Comment  |
|-------------|--------------------|--|
| (2000,0010) | Number of copies   | Any integer between 1 and 99. Default is 1                             |
| (2000,0020) | Print Priority     | Enumerated value: LOW, MED or HIGH. Default is MED. Attribute ignored. |
| (2000,0030) | Medium Type        | Enumerated value: PAPER, BLUE FILM or CLEAR FILM. Default is PAPER.    |
| (2000,0040) | Film Destination   | Enumerated value: MAGAZINE or PROCESSOR. Attribute ignored.            |
| (2000,0050) | Film Session Label | Up-to 64 characters.   |

The SCP will return one of the following status codes for N-CREATE:

| Code   | Status  | Meaning                              |
|--------|---------|--------------------------------------|
| 0x0000 | Success | Film session is successfully created |
| 0x0106 | Failure | Invalid attribute value              |
| 0x0213 | Failure | Resource limitation                  |
| 0xB600 | Warning | Memory allocation error              |

Note: ACETIAM PRI has a limitation for the number of created Film Sessions, and thus avoids too many films to be printed simultaneously.

➔ N-SET is used to update an instance of the Basic Film Session SOP Class.

The SCP will return one of the following status codes for N-SET:

| Code   | Status  | Meaning                              |
|--------|---------|--------------------------------------|
| 0x0000 | Success | Film session is successfully updated |
| 0x0106 | Failure | Invalid attribute value              |
| 0xB600 | Warning | Memory allocation error              |

➔ N-DELETE is used to delete the complete Basic Film Session SOP Instance hierarchy.



The SCP will return one of the following status codes for N-DELETE:

| Code   | Status  | Meaning                                    |
|--------|---------|--|
| 0x0000 | Success | Film session has been successfully deleted |

➔ N-ACTION is used to print the film session an instance of the Basic Film Session SOP Class.

The SCP will return one of the following status codes for N- ACTION:

| Code   | Status  | Meaning  |
|--------|---------|--|
| 0x0000 | Success | Films in the film session are accepted for printing  |
| 0xC600 | Failure | Film Session SOP Instance hierarchy does not contain Film Box SOP Instances.               |
| 0xB602 | Warning | Film Session SOP Instance hierarchy does not contain Image Box SOP Instances (empty page). |
| 0xC604 | Failure | Image position collision: multiple images are assigned to a single image position.         |
| 0x0213 | Failure | Resources limitation   |

#### 4.3.1.5 SOP Specific Conformance for Basic Film Box SOP Class

The following DIMSE services are supported:

N-CREATE

N-SET

N-DELETE

N-ACTION

➔ N-CREATE is sent by the SCU AE to create a Basic Film Box, once a Film Session has been successfully created.

The attributes of N-CREATE are shown in the following table:

| Tag         | Name                                 | Comment   |
|-------------|--------------------------------------|---|
| (2010,0010) | Image Display Format                 | STANDARD\m,n<br>ROW\r1,r2,...<br>COL\c1,c2,...  |
| (2010,0040) | Film Orientation                     | Enumerated value: PORTRAIT or LANDSCAPE. Default is PORTRAIT.   |
| (2010,0050) | Film Size ID                         | Enumerated value:<br>8INX10IN<br>10INX12IN<br>10INX14IN<br>11INX14IN<br>14INX14IN<br>14INX17IN<br>24CMX24CM<br>24CMX30CM<br>A3 or A4<br><i>See note</i>   |
| (2010,0060) | Magnification Type                   | Enumerated value: NONE, BILINEAR. Default is BILINEAR. Any other value will be considered and treated as BILINEAR.  |
| (2010,0080) | Smoothing Type                       | Attribute ignored   |
| (2010,0100) | Border Density                       | Enumerated value: WHITE or BLACK. Attribute ignored. Always BLACK.  |
| (2010,0110) | Empty Image Density                  | Enumerated value: WHITE or BLACK, or a valid Optical Density. Default is BLACK.   |
| (2010,0120) | Minimum Density                      | Default is 20 or may be found in GSDF calibration file  |
| (2010,0130) | Maximum Density                      | Default is 300 or may be found in GSDF calibration file   |
| (2010,0140) | Trim                                 | Default is NO. Attribute ignored.   |
| (2020,0500) | Referenced Presentation LUT Sequence | Accepted if the Presentation LUT SOP class has been negotiated, and a Presentation LUT SOP Instance has been successfully created.  |
| (0008,1150) | > Referenced SOP Instance UID        | Required if sequence is present.  |
| (0008,1155) | > Referenced SOP Class UID           | Required if sequence is present.  |
| (2010,015E) | Illumination                         | Accepted if the Presentation LUT SOP class has been negotiated, and a Presentation LUT SOP Instance has been successfully created. Default is 150, but may be set in the GSDF configuration file associated to the printer. |
| (2010,0160) | Reflected Ambient light              | Accepted if the Presentation LUT SOP class has been negotiated, and a Presentation LUT SOP Instance has been successfully created. Default is 0.  |

**Note:** ACETIAM PRI performs a mapping between DICOM defined terms and common paper sizes. This mapping overrides the input enumerated values and is user defined in ACETIAM PRI configuration menu.

The SCP will return one of the following status codes for N-CREATE:

| Code   | Status  | Meaning                                  |
|--------|---------|--|
| 0x0000 | Success | Film Box is successfully created         |
| 0x0120 | Failure | Missing one or more mandatory attributes |
| 0x0106 | Failure | Invalid attribute value                  |
| 0x0117 | Failure | Invalid object instance                  |
| 0x0110 | Failure | Processing failure                       |
| 0x0213 | Failure | Resources limitation                     |

➔ N-SET DIMSE service will allow updating the following Film Box attributes:

| Tag         | Name                      |
|-------------|---------------------------|
| (2010,0060) | Magnification Type        |
| (2010,0080) | Smoothing Type            |
| (2010,0100) | Border Density            |
| (2010,0110) | Empty Image Density       |
| (2010,0120) | Minimum Density           |
| (2010,0130) | Maximum Density           |
| (2010,0140) | Trim                      |
| (2010,0150) | Configuration Information |

The SCP will return one of the following status codes for N-SET:

| Code   | Status  | Meaning                          |
|--------|---------|----------------------------------|
| 0x0000 | Success | Film Box is successfully created |
| 0x0106 | Failure | Invalid attribute value          |
| 0x0110 | Failure | Processing failure               |
| 0x0105 | Failure | No such attribute                |

➔ N-DELETE is used to delete the Basic Film Box.

The SCP will return one of the following status codes for N-DELETE:

| Code   | Status  | Meaning                                |
|--------|---------|--|
| 0x0000 | Success | Film Box has been successfully deleted |

➔ N-ACTION is used to print one or more copies of a single film of the Film Session.

The SCP will return one of the following status codes for N- ACTION:

| Code   | Status  | Meaning  |
|--------|---------|--|
| 0x0000 | Success | Film accepted for printing   |
| 0xB603 | Warning | Film Box SOP Instance hierarchy does not contain Image Box SOP Instances (empty page). |
| 0xC604 | Failure | Image position collision: multiple images are assigned to a single image position.     |
| 0x0213 | Failure | Resources limitation   |

#### 4.3.1.6 SOP Specific Conformance for Basic Grayscale Image Box SOP Class

The following DIMSE services are supported:

##### N-SET

➔ N-SET may be used to update an instance of the Basic Grayscale Image Box.

The attributes of N-SET are shown in the following table:

| Tag         | Name                           | Comment  |
|-------------|--------------------------------|--|
| (2020,0010) | Image Position                 |  |
| (2020,0020) | Polarity                       | Enumerated value: NORMAL or REVERSE. Attribute ignored. Always NORMAL.   |
| (2010,0060) | Magnification Type             | Enumerated value: NONE, BILINEAR. Default is BILINEAR. Any other value will be considered and treated as BILINEAR. |
| (2010,0080) | Smoothing Type                 | Attribute ignored.   |
| (2020,0030) | Requested Image Size           |  |
| (2020,0110) | Basic Grayscale Image Sequence |  |
| (0028,0002) | >Samples Per Pixel             | 1  |
| (0028,0004) | >Photometric Interpretation    | Enumerated value: MONOCHROME1 or MONOCHROME2. Mandatory. No default.   |
| (0028,0010) | >Rows                          | Mandatory. No default.   |
| (0028,0011) | >Columns                       | Mandatory. No default.   |
| (0028,0100) | >Bits Allocated                | Valid range is 8 to 16. Mandatory. No default.   |
| (0028,0101) | >Bits Stored                   | Valid range is 8 to 16. Mandatory. No default.   |
| (0028,0102) | >High Bit                      | Valid range is 8 to 16. Mandatory. No default. Should be (Bits Stored - 1)   |
| (0028,0103) | >Pixel Representation          | Mandatory. No default  |
| (0028,0034) | >Pixel Aspect Ratio            |  |
| (7FE0,0010) | >Pixel Data                    | Mandatory. No default  |

The SCP will return one of the following status codes for N-SET :

| Code   | Status  | Meaning                                  |
|--------|---------|--|
| 0x0000 | Success | Image successfully stored in Image Box   |
| 0x0106 | Failure | Invalid attribute value                  |
| 0x0110 | Failure | Processing failure                       |
| 0x0120 | Failure | Missing one or more mandatory attributes |

#### 4.3.1.7 SOP Specific Conformance for Basic Color Image Box SOP Class

The following DIMSE services are supported:

##### N-SET

➔ N-SET may be used to update an instance of the Basic Grayscale Image Box.

The attributes of N-SET are shown in the following table:

| Tag         | Name                           | Comment  |
|-------------|--------------------------------|--|
| (2020,0010) | Image Position                 |  |
| (2020,0020) | Polarity                       | Enumerated value: NORMAL or REVERSE. Attribute ignored. Always NORMAL. |
| (2010,0060) | Magnification Type             | Attribute ignored. Always BILINEAR assumed.                            |
| (2010,0080) | Smoothing Type                 | Attribute ignored.   |
| (2020,0030) | Requested Image Size           |  |
| (2020,0110) | Basic Grayscale Image Sequence |  |
| (0028,0002) | >Samples Per Pixel             | 3  |
| (0028,0004) | >Photometric Interpretation    | RGB. Mandatory.  |
| (0028,0006) | >Planar Configuration          | 0 or 1   |
| (0028,0010) | >Rows                          | Mandatory. No default.   |
| (0028,0011) | >Columns                       | Mandatory. No default.   |
| (0028,0100) | >Bits Allocated                | 8  |
| (0028,0101) | >Bits Stored                   | 8  |
| (0028,0102) | >High Bit                      | 7  |
| (0028,0103) | >Pixel Representation          |  |
| (0028,0034) | >Pixel Aspect Ratio            |  |
| (7FE0,0010) | >Pixel Data                    | Mandatory. No default  |

The SCP will return one of the following status codes for N-SET:

| Code   | Status  | Meaning                                  |
|--------|---------|--|
| 0x0000 | Success | Image successfully stored in Image Box   |
| 0x0106 | Failure | Invalid attribute value                  |
| 0x0110 | Failure | Processing failure                       |
| 0x0120 | Failure | Missing one or more mandatory attributes |

#### 4.3.1.8 SOP Specific Conformance for Presentation LUT SOP Class

The following DIMSE services are supported:

##### N-CREATE

➔ N-CREATE may be used to create an instance of the Presentation LUT SOP Class.

The attributes of N-CREATE are shown in the following table:

| Tag         | Name                      | Comment                                     |
|-------------|---------------------------|---|
| (2050,0010) | Presentation LUT Sequence |   |
| (0028,3002) | > LUT Descriptor          |   |
| (0028,3003) | > LUT Explanation         |   |
| (0028,3006) | > LUT Data                |   |
| (2050,0020) | Presentation LUT Shape    | Enumerated values IDENTITY, INVERSE, LIN OD |

The SCP will return one of the following status codes for N-CREATE:

| Code   | Status  | Meaning                                |
|--------|---------|--|
| 0x0000 | Success | Image successfully stored in Image Box |
| 0x0105 | Failure | Unknown or unsupported attribute       |
| 0x0110 | Failure | Processing failure                     |
| 0x0111 | Failure | Duplicate SOP Instance                 |

#### 4.3.1.9 SOP Specific Conformance for Storage SOP Classes

The following DIMSE services are supported:

##### C-STORE

➔ C-STORE requests are used by the SCU to store datasets to ACETIAM PRI.

#### 4.3.1.10 Presentation Context Acceptance Criterion

N/A

#### 4.3.1.11 Transfer Syntax Selection Policies

ACETIAM PRI will always prefer Little Endian Explicit VR transfer syntax when provided by SCU clients. However, it may be configured to NOT support Little Endian Explicit VR.

## 5. Communication Profiles

### 5.1 Supported Communications Stacks

ACETIAM PRI provides DICOM V3.0 TCP/IP Network Communication Support as defined in PS 3-8 of the DICOM Standard.

### 5.2 TCP/IP Stack

ACETIAM PRI inherits its TCP/IP stack from the Windows system upon which it executes.

### 5.3 Physical Media Support

ACETIAM PRI is indifferent to the physical medium over which TCP/IP executes; it inherits this from the system upon which it executes.

## 6. Extensions/Specialization/Privatization

No extensions defined.

## 7. Configuration

ACETIAM PRI configuration is included in the application user interface through the Configuration Dialog.

### 7.1 AE Title/Presentation Address Mapping

Not Applicable.

### 7.2 Configurable Parameters

#### 7.2.1 Standard Configuration

ACETIAM PRI configurable parameters may be defined in the Configuration Dialog Box of the user interface. They are the following:

- TCP/IP port : default is 3100
- AE Title : default is PRINT
- Support of Basic Grayscale Print Management (on/off)
- Support of Basic Color Print Management (on/off)
- Support of Presentation LUT (on/off)
- Support of Storage to Print facility (on/off)
- Log events: Log all events or not (verbose) in the log window of the application.
- Debug: creates a detailed log file in a *logfiles* subdirectory besides the application.

- **Film Size Matching:** This option lets you specify the matching rules between DICOM defined Film IDs and common paper sizes.
- **GSDF calibration:** As ACETIAM PRI may manage several printers, you may define for each printer a GSDF calibration file that specifies: the default Illumination, the max DDL value, for a range of DDL input values, ranging from 0 to the maximum value, the corresponding measured ODs.

ACETIAM PRI has a limitation for the number of created Film Sessions, and thus avoids too many films to be printed simultaneously. Maximum number of simultaneous created Film Sessions is 6.

ACETIAM PRI may be configured to handle printing requests from client that have a strange (even illegal) interpretation of the standard. Two common cases are the following:

Print SCU negotiates successfully Grayscale and Color printing. Films boxes are then created using the grayscale presentation context for example, and thus grayscale image box instances are created and returned to the SCU. However, some SCU ignore the SOP Class of the returned image boxes, and do not hesitate to fill them with color pixel data using an N-SET command message on images. For this, in the N-SET message, the SCU even change the SOP Class UID of image to Color ImageBox SOP Class.

To handle this frequent illegal SCU behavior, ACETIAM PRI introduces internal "Proxy image boxes", that accept to dynamically change image box SOP Classes. ACETIAM PRI configuration panel allows making use of these "proxy" images boxes. This is however not the default.

Some Print SCU implementation of some well-known vendors will fail if an N-CREATE-RSP message for a FilmBox creation contains the optional Referenced Film Session information. ACETIAM PRI configuration panel let user change default SCP behavior, which is to NOT to return the Referenced Film Session Information.

## ***7.2.2 Optional Configuration***

### **7.2.2.1 Called AE Title routing**

According to ACETIAM PRI versions, it may possible to use the Called AE Title and the medium type to select a destination printer if several printers are available (local or network).

This particular configuration does not alter the general behavior of the software, and is described in a dedicated separate document (Administrator's Guide).

### **7.2.2.2 Calling AE Title limitations**

ACETIAM PRI may be configured to accept print requests from only a restricted set of Calling AE Titles.

This particular configuration does not alter the general behavior of the software, and is described in a dedicated separate document (Administrator's Guide).

## **8. Support of Extended Character Sets**

ACETIAM PRI supports Extended Character Set "ISO\_IR 100" Latin alphabet N° 1, supplementary set.

---



**NEHS DIGITAL**

1 rue Augustine Variot  
92240 Malakoff  
FRANCE

Tel: +33 (0)2 99 14 33 82

**ACETIAM Corp.**

162 Great Rd, Rear  
Acton, MA 01720  
USA

Toll Free: +1 877-384-2662 (USA & Canada)  
Phone/Fax: +1 617-953-0298